

## A Question of Taste – Demonstrator Notes

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Introduction – outside lab	<ul style="list-style-type: none"> <li>• Lab rules – use as questioning opportunity to warm up group</li> </ul>		0-10 min
Introduction	<ul style="list-style-type: none"> <li>• What will we be doing today – get ideas from the group</li> <li>• Summary of the day</li> </ul>		10-15 min
Recapping DNA background knowledge	<ul style="list-style-type: none"> <li>• Finding out about how much group knows about DNA and genes</li> <li>• Ask questions -               <ul style="list-style-type: none"> <li>○ What is DNA?</li> <li>○ What is a gene?</li> <li>○ Where is our DNA?</li> <li>○ What are the units of DNA called?</li> </ul> </li> </ul>		15-20 min
Tasting the PTC	<p><b>Students taste PTC and control strips</b></p> <p>Discussion of results</p> <ul style="list-style-type: none"> <li>• Ask what the students could taste – is it sweet, sour, bitter etc?</li> <li>• Distinction between bitter and sour</li> <li>• Ask who rated taste as a 5</li> <li>• Ask who rated taste as a 0</li> </ul>	<ul style="list-style-type: none"> <li>• Make sure students have washed hands</li> <li>• Ensure each student has a control (corner cut off) and PTC strip</li> <li>• Hand out record sheets</li> <li>• Tidy away strips and paper towels</li> </ul>	20-25 min

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Finding out about PTC	<ul style="list-style-type: none"> <li>• Explain that the taste strips were impregnated with the chemical phenylthiocarbamide (PTC)</li> <li>• Not everyone can taste this chemical</li> <li>• PTC tasting has a genetic basis – ask students what we call the visible expression of our genes? – phenotype</li> <li>• We are now going to investigate our genetics – what do we call the makeup of our genes? – genotype</li> </ul>		25 – 30 min
Harvesting cells	<ul style="list-style-type: none"> <li>• Explain to students that we are going to collect a cell sample by swilling mouth with saline solution (suggest that students gently chew the inside of their cheeks to improve their sample – we don't want any blood) – ask students why we use saline</li> <li>• Advise students not look at each other during swilling</li> </ul> <p><b>Students carry out mouthwash</b></p>	<ul style="list-style-type: none"> <li>• Make sure each student has saline sample</li> <li>• Give each student a number</li> <li>• Make sure students label 1.5ml microfuge tube with their number</li> <li>• Demo and supervise pipetting 1ml of sample into 1.5ml microfuge tube, make sure students take sample from bottom of cup</li> <li>• Collect tubes for centrifuging</li> </ul>	30 – 35 min
Separating out the cells	<ul style="list-style-type: none"> <li>• Ask students why we centrifuge the samples</li> <li>• Explain what will happen when samples are returned or if they have small/ fluffy pellets</li> <li>• Demonstrate tipping off most of the sample – point out level on the screen</li> <li>• Explain how to re-suspend sample</li> <li>• Ask students to place samples in tube rack for later</li> </ul>	<ul style="list-style-type: none"> <li>• Place tubes in centrifuge – place in order to speed returning</li> <li>• Check pellets before returning tubes to students – explain what to do if pellets are small</li> <li>• Supervise tipping off the supernatant</li> <li>• Supervise vortex mixing to re-suspend pellet</li> </ul>	35 – 45 min

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Micropipettes	<ul style="list-style-type: none"> <li>• Introduce idea of the precision needed to work with tiny volumes – ask students if they know what a micro litre is</li> <li>• Explain how to use the micropipettes – things to highlight –               <ul style="list-style-type: none"> <li>○ They are delicate and expensive!</li> <li>○ How to hold them correctly</li> <li>○ There are two stopping points that have different uses</li> <li>○ How to alter volume</li> <li>○ How to put on and remove a tip</li> <li>○ Avoiding contamination</li> </ul> </li> </ul> <p><b>Students practice pipetting</b></p>	<ul style="list-style-type: none"> <li>• Supervise pipette practice</li> <li>• Encourage students to measure a variety of volumes</li> <li>• Make sure students don't go outside the ranges of their pipettes – explain why this is</li> </ul>	45–55 min
Extracting DNA	<ul style="list-style-type: none"> <li>• Explain that 1.5ml microfuge tubes containing chelex have been placed in the trays</li> <li>• Demonstrate how to take 30µl of the cell suspension and transfer this to the chelex</li> </ul> <p><b>Students transfer sample</b></p> <ul style="list-style-type: none"> <li>• Ask students to collect samples in the rack and pass to the presenter ready to go in the heat block</li> </ul>	<ul style="list-style-type: none"> <li>• Hand out chelex – 1 sample per student</li> <li>• Supervise transfer of cell suspension</li> <li>• Collect samples from group and place in 100°C heat block – set timer for 10 minutes</li> </ul>	55 min – 1 hr
While samples are in heat blocks	<ul style="list-style-type: none"> <li>• Explain why we use chelex – use <i>PowerPoint slide 20</i></li> <li>• Ask students why we are heating the samples to 100°C</li> <li>• Ask why we will centrifuge the samples again</li> </ul>	<ul style="list-style-type: none"> <li>• Remove samples from heat block after 10 minutes</li> <li>• Spin samples in centrifuge for 1 minute</li> </ul>	1hr – 1hr 10

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Taste – how it works and the <i>TAS2R38</i> gene (during 10 minute heat block)	<ul style="list-style-type: none"> <li>Explain that taste receptors are proteins</li> <li>PTC taste receptor coded for by <i>TAS2R28</i> on chromosome 7</li> </ul>	<ul style="list-style-type: none"> <li>During this time samples will finish in heat block/ centrifuge so return samples to students</li> </ul>	
PCR set up	<ul style="list-style-type: none"> <li>Ask students if they remember what the process we are doing today is called</li> <li>Explain that we will think about how it all works once we have set up the reaction</li> <li>Explain the different ingredients needed – see if the students can name any of them</li> <li>Ask students to label PCR tubes with numbers</li> <li>Demonstrate how to take DNA sample - <i>use PowerPoint slide 26 to show where to concentrate</i></li> <li>Explain how much to add – <i>see PowerPoint slide 27</i></li> </ul>	<ul style="list-style-type: none"> <li>Hand out tubes with primers (one per pair) and PCR beads (one each)</li> <li>Supervise students taking their DNA samples</li> <li>Make sure students take samples from the top of the liquid avoiding beads etc</li> <li>Supervise mixing of ingredients</li> <li>Collect PCR tubes and place in PCR machine</li> </ul>	1 hr 10 – 1hr 20
PCR background	<ul style="list-style-type: none"> <li>Ask students if they can explain what PCR does</li> <li>Explain that it can be used to copy specific areas of DNA</li> <li>Show flash animation – describe each step</li> </ul>		1hr 20 – 1hr 35
PCR activity	<ul style="list-style-type: none"> <li>Explain how the activity works</li> </ul> <p><b>Students carry out activity</b></p>	<ul style="list-style-type: none"> <li>Hand out the following to groups of 4 students <ul style="list-style-type: none"> <li>1 piece of double stranded DNA</li> <li>Some primers</li> <li>Some nucleotides</li> <li>Scissors</li> <li>tape</li> </ul> </li> <li>Supervise students as they carry out the activity</li> </ul>	1hr 35 – 1hr 50

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Gel pouring	<ul style="list-style-type: none"> <li>Briefly explain that later we will be thinking about electrophoresis but now we need to start setting it up</li> <li>Ask for volunteers from each table</li> <li>Demonstrate how to place combs, casting gates and pour gel</li> </ul> <p><b>Students help prepare gels</b></p>	<ul style="list-style-type: none"> <li>Hand out gels to students – make sure they take care of hot gel</li> <li>Supervise gel pouring</li> </ul>	1hr 50 – 2hr
<b>Lunch</b>			2hr – 2hr 30
How to work out if you have taster DNA	<ul style="list-style-type: none"> <li>Introduce the idea of SNPs – try to get students to work out what a SNP might be if they don't already know</li> <li>Ask students what we call different versions of the same gene</li> <li>Explain the taster and non-taster SNP</li> <li>Explain how a single base change may or may not affect the protein coded for</li> <li>Try to get students to work out why there seems to be a whole range of tasting abilities not just taster and non-taster (<i>explored in PowerPoint slide 38 if students can't answer</i>)</li> </ul>		2hr 30 – 2hr 50
Homozygous and heterozygous	<ul style="list-style-type: none"> <li>Explain that we have two copies of <i>TAS2R38</i> and that the various combinations of the C145G SNP explains the range of tasting ability</li> <li>Explain the terms homozygous and heterozygous</li> </ul>		

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Explaining restriction enzymes	<ul style="list-style-type: none"> <li>• Ask if the students know what restriction enzymes are</li> <li>• Explain how the enzymes cut DNA at specific places</li> <li>• Explain HaeIII and its restriction site</li> <li>• Show students the two versions of the PTC receptor gene (<i>slide 42</i>) – can they work out which one will be cut by the HaeIII?</li> </ul>	<ul style="list-style-type: none"> <li>• Get ice beakers from the freezer</li> </ul>	
Setting up the digest	<ul style="list-style-type: none"> <li>• Ask students to label fresh tube with number and the letter D</li> <li>• Explain how to set up the digest – see <i>PowerPoint slide 43</i></li> <li>• Explain that we will set up a negative control by heating the remaining PCR product with no HaeIII</li> </ul> <p><b>Students set up restriction digest</b></p>	<ul style="list-style-type: none"> <li>• Hand out tubes of HaeIII, water</li> <li>• Collect all tubes and place in 37 degree heat bath/ block for 30 min</li> </ul>	2hr 50 – 3hr
Projected digest results	<ul style="list-style-type: none"> <li>• Run through the possible digest results - <i>PowerPoint slide 44 (note – the diagram in this slide represents our PTC product, not our whole chromosome)</i></li> </ul>		

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Restriction enzyme/ electrophoresis activity	<ul style="list-style-type: none"> <li>• Explain how the activity works (see instructions)</li> </ul> <p>Things for students to know –</p> <ul style="list-style-type: none"> <li>• The double stranded DNA will have the same restriction enzyme site with the same sequence which reads forwards on the sense strand and backwards on the antisense</li> <li>• Restriction enzymes only cut at specific sites</li> <li>• Both strands are cut</li> <li>• We are using just one of many restriction enzymes</li> <li>• Other enzymes will cut different sites</li> </ul> <p><b>Student carry out restriction enzyme activity</b></p> <p>Students can now simulate electrophoresis with their DNA fragments – either on the table top or wall etc</p> <p>Things for students to know –</p> <ul style="list-style-type: none"> <li>• Small fragments move further than large fragments</li> <li>• Different DNA samples will produce different banding patterns</li> <li>• This technique can be used for DNA fingerprinting etc</li> </ul>	<ul style="list-style-type: none"> <li>• Hand out the following <ul style="list-style-type: none"> <li>○ One DNA strand per group</li> <li>○ Scissors</li> <li>○ Velcro, blue tack etc if sticking on the wall or board</li> </ul> </li> </ul>	3hr – 3hr 20

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Humans, chimps and PTC	<ul style="list-style-type: none"> <li>• Remind students that in human populations there are PTC tasters and non-tasters</li> <li>• Ask them if they think chimps can taste PTC, why do they think this?</li> <li>• Ask why humans and chimps may have evolved similar traits – useful to both species</li> <li>• Explain the background research looking at evolution of non-tasting in humans and chimps</li> <li>• Show the different SNPs and explain that these haven't been inherited from the same place, but the outcome is the same</li> <li>• What do we call this? Convergent evolution <ul style="list-style-type: none"> <li>○ Get students to come up with other examples of convergence on their tables</li> </ul> </li> </ul>		3hr 20 – 3hr 35
Gel electrophoresis	<ul style="list-style-type: none"> <li>• Briefly recap how electrophoresis works</li> <li>• Ask for volunteers from each table to help prepare the gels – removing combs, casting gates and pouring the buffer</li> </ul>	<ul style="list-style-type: none"> <li>• Supervise gel preparation, making sure students are careful and don't damage gels</li> <li>• Hand out bottles of buffer</li> <li>• Collect up combs, casting gates</li> </ul>	3hr 35 – 3hr 45
Loading samples	<ul style="list-style-type: none"> <li>• Demonstrate how to load the samples</li> <li>• Explain that everyone will have a chance to practice loading</li> <li>• Explain that we will load both our digested and undigested samples – show lane sheet for recording</li> </ul> <p><b>Students load their practice and real samples</b></p>	<ul style="list-style-type: none"> <li>• Return samples from heat block</li> <li>• Supervise loading – spin samples if needed</li> <li>• Set up and turn on power packs (120V for 25 min)</li> </ul>	3hr 45 – 4hr 10

Section	Learning outcome/ things that need to be covered	Demonstrator/ tech activity	Timing
Dialogue	<p>Explain the dialogue activity –</p> <ul style="list-style-type: none"> <li>The activity is designed to look at the possible reasons for why non tasters might have evolved in human and chimp populations</li> <li>No right or wrong answer</li> <li>Non-tasting may be a coincidence but there is more likely a reason – i.e. an advantage</li> </ul> <p>Students will –</p> <ul style="list-style-type: none"> <li>Arrange prompt cards in order from good reason to bad (cards can go side by side if they want)</li> <li>Have to be able to explain their reasoning</li> </ul> <p>Activity works well if started in small groups before then having a class discussion</p>	<ul style="list-style-type: none"> <li>Hand out student cards</li> <li>General tidy of lab space</li> <li>Put out light boxes and filter glasses</li> </ul> <ul style="list-style-type: none"> <li>Towards the end of the dialogue the electrophoresis will finish – remove gels and place on light boxes</li> </ul>	4hr 10 – 4hr 40
Results	<ul style="list-style-type: none"> <li>Explain that the gels have now finished</li> <li>Describe how the stain works</li> <li>Run through the possible results we will see and what these tell us</li> <li>Allow students to look at their gels</li> </ul>	<ul style="list-style-type: none"> <li>Help students interpret their gels</li> </ul>	4hr 40 – 4hr 55
Round up	<ul style="list-style-type: none"> <li>Ask students if their genotype matched their phenotype</li> <li>If any don't match try and get students to explain why this might happen</li> </ul>		
End of workshop	<ul style="list-style-type: none"> <li>Ask students/ teachers to fill in evaluation forms</li> <li>Thank students for their work</li> <li>Ask students to wash hands and return lab coats</li> </ul>		4hr 55 – 5hr